97050748  
rev. 002  
07/2014

**iM3<sup>®</sup>**

Revolution 4DC

EN  
IT  
FR  
DE  
中文

# OPERATING INSTRUCTIONS

## CONTENTS

1. General precautions .....	4
1.1. Description of the symbols used .....	4
1.2. Purpose and use of the equipment .....	5
1.2.1. Classification .....	5
1.2.2. Environmental conditions .....	5
1.2.3. Guarantee.....	5
2. Description of the system .....	6
2.1. Description of the x-ray unit .....	6
2.2. Operating mode .....	6
3. Operation.....	7
3.1. Turning the x-ray unit on and off .....	7
3.1.1. Turning on the basic X-ray unit.....	7
3.1.2. Turning on the handheld.....	7
3.1.3. Automatic handheld shut off .....	8
3.1.4. Handheld stand-by according to time .....	8
3.2. Handheld display functions .....	8
3.3. Control panel.....	8
3.4. Checking the parameters.....	9
3.5. Factory settings.....	10
3.6. Batteries and charge level indication .....	10
3.7. X-ray generator indicator light.....	10
4. Using the x-ray unit .....	11
4.1 Position of the patient .....	11
4.2. Positioning the x-ray head .....	11
4.2.1. Hypersphere technology.....	11
4.3. Setting the exposure mode and time .....	12
4.4. Setting the mode and exposure time in USER mode .....	12
4.5. Procedure to be followed when taking the x-ray .....	13
5. Advanced options .....	14
5.1. Setting the safety unlock mode .....	14
5.2. Setting the operating mode .....	14
5.3. Restoring factory settings .....	15
6. Error messages .....	15
7. Routine maintenance.....	16
8. Cleaning and disinfection .....	16
9. Dismantling the equipment when no longer used.....	17
10. Technical specifications .....	17
11. Dimensions .....	18
12. Identification nameplates.....	20

## 1. General precautions

These instructions explain how to properly use the x-ray unit. Please carefully read and become familiar with the content of this manual before attempting to use the apparatus.

 **NOTE:** This manual does not specify the rules and regulations for acquiring, possessing or using a source of ionizing radiation as each country has its own laws. Only the most common ones shall be mentioned and this means that it is the user's responsibility to check local standards and observe the relevant laws.

No part of this manual is to be reproduced, stored in a retrieval system or transmitted in any form or by any means, i.e. electronic, mechanical, photocopying, translation or otherwise, without the prior written permission of the manufacturer.

CEFLA SC - CEFLA DENTAL GROUP has a policy of ongoing development and is committed to improving its products. Therefore, some instructions, specifications and illustrations given in this manual may differ slightly from the product actually purchased. The manufacturer reserves the right to make changes to this manual without giving prior notice.

The original version of this manual is written in Italian. Translation from the original in Italian

### 1.1. Description of the symbols used

The following symbols are used:

	<b>WARNING!</b> Indicates a situation where failure to observe the instructions could cause damage to the equipment or injury to the user and/or patient.
	<b>NOTE:</b> Indicates important information for the user and/or technical support personnel.
	Earth connection.
	Alternating current.
	On.
	Off.
	ionizing radiation.
	Equipment abides by the essential requirements established by USA and Canada.
FCC ID	F.C.C. mark (Federal Communication Commission)

## 1.2. Purpose and use of the equipment

This x-ray unit is designed for use in the dental veterinary surgery to make endo-oral x-rays for diagnostic purposes. This equipment can be used to produce traditional x-rays developed using chemicals or, alternatively, it can be used with digital x-ray sensors.

### **ATTENTION!**

This x-ray unit is **FOR VETERINARY USE ONLY**.  
This x-ray unit is **NOT FOR HUMAN USE**.

### 1.2.1. Classification

**RADIO EQUIPMENT AND TELECOMMUNICATIONS TERMINAL EQUIPMENT** Classification  
Equipment classification in accordance with directive 99/05/EEC art.12: Class I.

#### **EMC classification**

Classification of the equipment according to standard CEI EN 55011: **Group 1 - Type B**.

### 1.2.2. Environmental conditions

The equipment should be installed in an area where the following conditions are present:

- Temperature from +10° to +40° C;
- relative humidity from 25 to 75% without condensate;
- atmospheric pressure from 700 to 1060 hPa;
- The electrical wiring in the room where the equipment is installed must conform to the I.E.C. 60364-7-710;V2 specification (i.e. the regulations concerning the electrical wiring to be used in surgeries) or equivalent standards in force in the country where the equipment is installed.
- **ELECTRICAL CONNECTIONS** : The electrical wiring must have an effective ground conductor as set forth by I.E.C. – US National Electrical Codes and C.E.I. standards. In Italy, electrical wiring must comply to standards IEC 60364-7-710 which require that a ground fault circuit interrupter is installed upstream. The ground fault circuit interrupter must be as specified below:
  - contact capacity: 250 V 10A in compliance with standards IEC 60898-1 and IEC 60947-2.
  - differential sensitivity: 0.03A.
  - power supply: 3 x 2.5 sq.mm

The color of the 3 wires should be as specified in the standards (brown power, BLUE neutral, YELLOW/GREEN ground).

### 1.2.3. Guarantee

CEFLA SC - CEFLA DENTAL GROUP guarantees the safety of the operator as well as the reliability and high performance of the equipment. The guarantee shall be voided if the following precautions are not taken:

- Observe the conditions specified in the guarantee certificate itself.
- The equipment must only be used by following the instructions contained herein.
- Equipment installation, expansion and technical support must be performed exclusively by personnel authorized by the manufacturer to carry out these operations.
- Never open the equipment casing. Installation, repairs and, in general, any other operations requiring the casing to be opened are to be performed exclusively by personnel authorized by the manufacturer to carry out these operations.
- The equipment must be installed in areas which comply with that stated above at point 1.2.2 “Ambient Conditions”.
- The area where the x-ray unit is installed must comply with official regulations regarding protection against radiation in the country where the equipment is used.

### **SAFETY PRECAUTIONS**

- If any person who is not an authorized technician changes the product in any way by replacing parts or components with other ones not used by the manufacturer they shall assume responsibility for the product.
- Do not forget to turn off the main switch on the equipment before leaving the surgery.
- This equipment is not water-proof (risk of electrocution).
- This equipment is not suitable for use in the presence of a mix of inflammable anaesthetic gas with oxygen or nitrous oxide.
- This equipment must be stored properly so that it is kept in top working order at all times.
- Use of electric scapels or other electric apparatus that do not comply to standard I.E.C. 60601-1-2, in the doctor's office or in the immediate vicinity of it may cause electromagnetic interferences or interferences of other types, resulting in equipment malfunctions. In these cases shut off the power supply to the equipment before hand.
- The manufacturer shall not be held responsible for misuse, carelessness or improper use of the equipment.
- This equipment is to be used exclusively by qualified personnel with the proper medical training.
- The user must be present at all times when the equipment is turned on or ready for start-up. In particular, never leave the equipment unattended in the presence of children/the mentally disabled or other unauthorised personnel in general.
- If the x-ray equipment is damaged or oil leaks, do not use the equipment and contact customer service immediately.

### **Protection against radiation**

**X-rays are hazardous and adequate precautions must be taken when using them. Areas where it is possible to be exposed to the x-rays shall be clearly indicated by using this symbol which should remind personnel to observe the safety rules laid down by the laws in force in the country where the equipment is used.**

- Control the emission of x-rays from the greatest distance possible (several metres) from the focal spot and the X-ray irradiation beam in the opposite direction to where the rays are emitted.
- Only the authorized personnel and the patient can remain in the area when x-rays are being emitted.
- Always protect the patient's thyroid and gonads under all circumstances.

## 2. Description of the system

### 2.1. Description of the x-ray unit

#### Description of the parts:

##### a X-ray generator.

Depending on the operating mode, the high frequency constant potential x-ray generator works at 60KV 7ma (En60 mode, 63KV 6ma (En63 mode) or 65KV 6ma (En65 mode).

The generator can turn on a horizontal plane endlessly; on the other hand, as far as vertical movement is concerned, rotation is limited upwards by a mechanical end-stop.

##### b Removable cone.

The generator can work with different types of collimators that are automatically recognized:

- 8" ROUND COLLIMATOR (incorporated in the generator): minimum skin/focus distance 20 cm and 60 mm output beam
- REMOVABLE 12" round CONE: minimum source/skin distance 30cm and diameter of x-ray beam exiting cone 55mm (with cone attached).

##### c Focus spot.

##### d Double pantograph arm.

##### e Extension arm.

The extension is available in three lengths: 40 cm, 60 cm and 90 cm.


##### f Handheld.

The handheld can be placed either near the control unit or in a remote position. As a result, the doctor can move conveniently around the room and move out of the area where x-rays are emitted.

##### g Hand held holder.

##### h Control unit.

##### i Master switch.

### 2.2. Operating mode

The x-ray unit can function in different modes (MULTI-MODE technology):

#### 1) En60 operating mode

X-ray emission at 60KV and 7ma. The X-ray unit automatically calculates the best exposure time (from 0.01s to 1.00s) based on the selected tooth and patient size.

#### 2) En63 operating mode

X-ray emission at 63KV and 6ma. The X-ray unit automatically calculates the best exposure time (from 0.01s to 1.00s) based on the selected tooth and patient size.

#### 3) En65 operating mode

X-ray emission at 65KV and 6ma. The x-ray unit automatically calculates the best exposure time (from 0.01s to 1.00s) based on the selected tooth and patient size.

#### 4) AUTO operating mode

The x-ray unit automatically suggests the best operating mode (En60, En63 or En65) and exposure time (from 0.01s to 1.00s) based on the selected tooth and patient size.

The suggested exposure time can be corrected from the handheld.

In addition, the x-ray unit has a special USER mode that can be set by the user. The user can select the best combination of load factors (operating mode and exposure time) for each tooth and patient size.

### 3. Operation

#### 3.1. Turning the x-ray unit on and off


##### 3.1.1. Turning on the basic X-ray unit

The control unit is turned on and shut off with switch (A).  
The switch lights up to signal the control unit is energized.

 **NOTE:** The technical specifications of the switch are given in paragraph 1.2.2.

The control unit is turned on and off with the main switch (A), as illustrated in the figure below. The switch lights up when the control unit is energized. Whenever turned on, the equipment performs an operational test that takes a few seconds. A beep is provided at the end of the test.

 **NOTE:** The exposure time and the parameters displayed when the unit is turned on are the last ones set before the central control unit was turned off. If the central control unit is left untouched for a few minutes it will go into stand-by mode. Simply press any key on the control panel to reactivate it.

##### 3.1.2. Turning on the handheld

The handheld is turned on by pressing any key, except for the one for x-ray emission. A buzzer rings to confirm the apparatus has been turned on. The unit will be in the standard configuration described in detail in paragraph 3.1.3 and then search for the base it works with.

If the base is off, the handheld will not indicate the field nor status "ready". If the base is latter turned on, the handheld will detect it within thirty seconds or by pressing any function key on the pushbutton panel.

 **NOTE:**  
To optimize the range of the handheld while it is being used, keep it away from walls and metal instruments and above all, do not cover its antenna found on top of the screen. In addition, performance may be reduced if the handheld is moved too quickly while x-rays are being taken. Error E 31 may be displayed if out of range.


### 3.1.3. Automatic handheld shut off

Once the control unit has been turned off the handheld automatically shuts off after approximately one minute. The handheld also automatically shuts off when it is further than the maximum allowable distance away from the control unit.

### 3.1.4. Handheld stand-by according to time


The entire x-ray unit will switch over to stand-by (even if the base is on) and the handheld will automatically shut off after approximately five minutes of non-use to save battery power. The handheld turns back on displaying the last selection made by the user whenever any key, except for the X-ray emission key, is pressed. To edit the stand-by time, refer to chapter 4 that deals with the handheld's "Advanced options".

## 3.2. Handheld display functions

## 3.3. Control panel

As illustrated in the figure below, the handheld has four function keys and a single x-ray emission key.

The main functions of the keys on the handheld vary according to how they are pressed:

KEY	BRIEFLY PRESSED (less than 3 sec.)	PRESSED LONGER (more than 3 sec.)
	Changes over from LARGE to SMALL and vice versa (takes place when key is released).	Saves the selected setting (exposure time, sensitivity, etc...). The memo icon () lights up when the data item can be saved.
	Changes amongst the various types of teeth to select the area to be examined.	It displays the values corresponding to the tooth exposure times in mGy and in mGy*cm2 if pressed again.
	Increases the exposure times in steps, according to the set scale.	Increases the scroll speed of the values in increasing order.
	Increases the exposure times in steps, according to the set scale.	Increases the scroll speed of the values in decreasing order.
	NO EFFECTS ARE OBTAINED IF THE KEY IS PRESSED LESS THAN A SECOND.	STARTS X-RAY EXPOSURE (the button has to be held down while the x-rays are being emitted, "dead man" function).


- NOTE:** "Dead man" function: the system that starts x-ray exposure with the dedicated key on the wireless handheld allows x-rays to be emitted only when the user presses and holds down the exposure key. X-ray emission will stop if the key is released ahead of time.
- NOTE:** The function related to pressing the key briefly is performed by pressing the key which will activate the function assigned to it. On the other hand, to perform the function carried out when the key is held down longer, press the key until the relative function is started. The buzzer will ring shortly to signal the function has started.
- NOTE:** Warm-up: When the equipment has not been used for a prolonged period (more than 3 months) or when turned on for the first time, a number of emissions with short times (0.01-0.02 sec.) are recommended and then some pictures with 0.1 sec. intervals to better stabilize operation of the x-ray tube before using it.

### 3.4. Checking the parameters

Before actually taking an exposure, make sure the exposure parameters for the examination in progress are correctly set.

- **Checking the selected body build.**
  - "SMALL" selected: indicates the x-ray unit is set for patients with small builds.
  - "LARGE" selected: Indicates the x-ray unit is set for patients with average-large builds.

**NOTE:** After the change has been made, the preset exposure times will automatically be modified.

• **Checking the selected type of intraoral exam.**

 <p>Upper molars</p>	 <p>Lower incisors</p>
 <p>Upper canines/bicuspids or rear "bitewing"</p>	 <p>Lower canines/bicuspids</p>
 <p>Upper incisors or front "bitewing"</p>	 <p>Lower molars</p>

### 3.5. Factory settings

The x-ray unit is supplied with the following factory settings:

- Operating mode: AUTO.
- Sensitivity: level 19.
- Handheld stand by: 5 minutes
- Exposure times as per standard R20: 0,010 - 0,011 - 0,012 - 0,014 - 0,016 - 0,018 - 0,020 - 0,022 - 0,025 - 0,028 - 0,032 - 0,036 - 0,040 - 0,045 - 0,050 - 0,056 - 0,063 - 0,071 - 0,080 - 0,090 - 0,100 - 0,110 - 0,125 - 0,140 - 0,160 - 0,180 - 0,200 - 0,220 - 0,250 - 0,280 - 0,320 - 0,360 - 0,400 - 0,500 - 0,560 - 0,630 - 0,710 - 0,800 - 0,900 - 1,000.

 **NOTE:** These times comply with current standards I.E.C. 60601-2-7 (1999) and the ISO 497 series R'20 recommendations and **CANNOT BE MODIFIED.**

### 3.6. Batteries and charge level indication

The handheld runs on two widely available AA alkaline batteries to assure sufficient stand-alone operation. The charge level of the batteries is given on the screen as follows:


-  Battery fully charged (a symbol does not appear in the area that shows the battery charge level).
-  Battery half-charged.
-  Battery charge level low or almost dead (causing the handheld to automatically shut off).

 **NOTE:** The batteries should be removed from the handheld if it is not going to be used for an extended period.

### 3.7. X-ray generator indicator light

The x-ray generator comes with an indicator light (B) that signals apparatus status:

- color (purple) > x-ray unit on (regular condition)
- Flashing color (purple) > stand-by (low consumption)
- Blue > x-ray one – head released
- Yellow > x-rays being emitted
- Red > fault

## 4. Using the x-ray unit

### 4.1 Position of the patient

A positioner or alignment device specific for the selected image receiver should always be used to assure the x-rays are correctly aligned regardless of the position the patient's head is in.

### 4.2. Positioning the x-ray head


Position the x-ray head so that the cone is aligned with the image receiver.

#### 4.2.1. Hypersphere technology

Hypersphere technology allows the x-ray head to turn endlessly on both the horizontal and vertical planes.

The x-ray head is initially blocked by an electromechanical brake.

The head can be tilted into the position required to take the x-ray by touching the unlock areas. To lock it again, release the unlock areas.

 **NOTE: Firmly hold the head with both hands when putting it in place.**

It is possible to set a safety unlocking mode that allows the head to be turned only by pressing both unlock buttons. This prevents the head from unlocking unexpectedly after one of the two unlock buttons has been accidentally pressed. To activate this mode, refer to "Advanced options" in chapter 5.

## 4.3. Setting the exposure mode and time

The exposure parameters are set by following the directions given below:

- 1) select the tooth to be examined
- 2) select the patient size

The exposure time is automatically shown on the handheld screen.

 **NOTE: each tooth and patient size selected is displayed for approximately 1 second according to the operating mode (En60, En63 or En65) used.**

The suggested exposure time can be changed with keys  and . Exposure times ranging from 0.01s and 1.00s belonging to the R'20 scale can be set. Random exposure times different from the ones provided in the R'20 scale cannot be set.

When the exposure time displayed differs from the default setting, icon  comes on.

To save the new setting, make sure icon  is on and then press and hold down key  for approximately 2 seconds. The handheld will beep to confirm the setting has been saved. At this point, make sure icon  is off.

 **NOTE: if the exposure time is not saved, the change made will be lost after a new entry or as soon as the handheld changes over to stand-by.**

 **Important: after customized settings have been made, the “Original exposure values charts” are no longer valid.**

If icon  is displayed while the exposure time is changed, it means the set time cannot be saved for the selected tooth-patient size combination. In any case, the x-rays can be taken with the set time.

 **Important: when the suggested exposure time is changed, the sensitivity factor is also modified (by default set to F=19). Once this change has been saved, it is applied to all the teeth and both patient sizes**

The exposure time can also be modified by changing the sensitivity factor. Press keys  and  at the same time, the actual sensitivity factor will be displayed.

Use keys  and  to change the value from 3 to 25. If the displayed value differs from the one previously saved, icon  comes on. To

quit this mode, press key  or . The change made to the sensitivity factor is applied to all the teeth and both patient sizes.

The selected operating mode is always used for each tooth and patient size combination in modes En60, En63 and En65.

In AUTO mode, each tooth and patient size combination is associated to the best mode from amongst the ones available. In this mode it is not possible to assign a mode other than the default one to each combination. To set the mode, refer to paragraph 4.5 “Setting the mode and exposure time in USER mode”.

To change the mode amongst En60, En63, En65 and AUTO refer to paragraph 5.2 “Setting the operating mode”.

## 4.4. Setting the mode and exposure time in USER mode

In USER mode, it is possible to assign an exposure time and a mode from amongst En60, En63 and En65 to each tooth-patient size combination. The default setting corresponds to the AUTOMode settings with sensitivity factor F=19

To activate USER mode regardless of the mode currently being used, press keys  and  at the same time. Icon  will come on to signal USER mode is active.

To deactivate USER mode press keys  and  again (icon  goes off).

The exposure parameters are set as directed below:

- 1) select the tooth under examination
- 2) select the patient size.

The exposure time is automatically displayed on the handheld.

 **NOTE: It is not possible to access the sensitivity factor menu in USER mode. In addition, keys  and  are inoperative in User mode.**

The exposure times and mode assigned to the tooth – patient size combinations are custom set by following the directions given below:

- 1) press and hold down key  about two seconds. Customized settings can be entered and icon  comes on.
- 2) select the desired tooth-patient size combination
- 3) change the exposure time with keys  and .

 **NOTE: it is possible to set exposure times ranging from 0.01s and 1.00s that are part of the R'20 scale.**

- 4) press keys  and  simultaneously to open the menu used to select the operating mode
- 5) select the operating mode with keys  and .
- 6) quit the menu and press key  to make the entry operative (if key  is pressed, the menu will be quit without changing the previous setting).
- 7) press and hold down key  for approximately two seconds to confirm the entry and disable customized settings (icon  goes out).

 **NOTE: it is possible to set the exposure parameters for several combinations. To do this, repeat steps 2 to 6 before going on to step 7.**

## 4.5. Procedure to be followed when taking the x-ray

- Pick up the handheld and go a safe distance away (at least 2 meters) maintaining visual contact with the patient and x-ray unit during the exposure. Make sure “READY” is indicated.

# READY

- Tell the patient to stay still.
- Press and hold down the “Exposure” key on the handheld until the audible warning sound (beep) stops and the yellow light goes out.


“X-ray emission” key

Light on control panel illuminated during x-ray emission

 **NOTE: If the “EMIT X-RAY” key is released at any time, exposure will be interrupted and error code E01 will appear on the display.**

- Once exposure has been completed, it is possible to proceed with the next exposure unless the x-ray unit has reached the maximum allowable temperature. The percentage the cone exceeds the maximum allowable temperature is always shown on the screen (see icon below).

- Once the temperature has been reached, wait the pause time for cooling signaled by symbol .
- At this point the exposure function will be disabled until the screen shows "READY" again.
- As soon as "READY" appears on the handheld, another exposure can be taken.

## 5. Advanced options

The handheld allows the user to view, edit and set some operating parameters by simply combining the keys provided. Follow the steps given below to access:

Key combination	Description of command
 + 	Press these two keys to adjust the sensitivity levels (determined based on the table given below and type of sensor/receiver used), modifying the current value from the minimum to the maximum allowable one (on a scale from 3 to 25), with keys "+" and "-". Press key "size" to confirm the desired level and go back to the main screen. This menu is not available in USER mode.
 + 	Hold down these two keys to go to the set up menu (from P 01 to P 07). Press key "Build" to make the selection. Once within the individual configurations, they can be scrolled with keys "+" and "-" and selected by pressing key "Build" again. Key "tooth" quits set up without saving the setting. The configurations are given in detail below: P 01: Sets the stand by time (from a minimum of 5 to a maximum of 30 minutes). P 02: Assigns an identification tag to the x-ray unit's base (from 1 to 5 or none). P 03: Shows the list of software versions- P 04: Handheld code display P 05: Activates/deactivates the safety unlock mode (see section 5.1) P 06: Selects the operating mode (En60, En63, En65 and AUTO). P 07: Sets the type of removable cone used
 + 	Activating/deactivating the USER mode Icon  comes on to signal USER mode is activated.

### 5.1. Setting the safety unlock mode

The x-ray unit has a safety unlock for the ball joint.

The default setting allows the ball joint to be disengaged by simply touching one of the keys present on the front of the head. To prevent accidental contact with the keys from unexpectedly disengaging the ball joint (and therefore causing undesired movement of the head), the safety unlock mode can be activated. In this mode, the ball joint is disengaged only if both keys are activated at the same time.

To set the safety unlock mode, press keys  and  to go to the set up menu. Scroll the parameters up to parameter P05 and press key . Scroll the options to select "ON" and press key . Press key  to quit the set up menu.

### 5.2. Setting the operating mode

The x-ray unit features the following operating modes:

- **En60:** all the x-rays are taken at 60KV and 7mA.
- **En63:** all the x-rays are taken at 63KV and 6mA.
- **En65:** all the x-rays are taken at 65KV and 6mA.
- **AUTO:** the system automatically selects the best setting from amongst En60, En63 and En65 for each tooth-patient size combination.

 **NOTE: the current setting is displayed on the handheld for approximately 1 second for each tooth-patient size selected before the relative exposure time is shown.**

To set the operating mode, press keys  and  to go to the set up menu. Scroll the parameters up to parameter P06 and then press key . Scroll the options to find the desired operating mode and then press key . Press key  to quit the set up menu.

### 5.3. Restoring factory settings

To restore the factory settings (see paragraph 3.5) press keys  and  to go to the set up menu. Press keys  and  simultaneously. "rESS" will briefly appear and the handheld will be rebooted.

### 6. Error messages

 **NOTE:** In presence of strong wireless communication traffic, the connection between handheld and generator may break down. To restore the connection run the "Restoring factory settings" procedure.

FAULT	CAUSE	SOLUTION
<b>E01</b>	X-RAY KEY RELEASED TOO EARLY	Hold down the key until the image has been captured.
<b>E02</b>	SHOOTING SEQUENCE NOT COMPLETED	Handheld most likely lost the signal. Try to repeat exposure. If the problem persists, contact technical service.
<b>E03</b>	HANDHELD INTERNAL TEST ERROR	Take out the batteries and then put them back in after waiting a few seconds. If the problem persists, contact technical service.
<b>E04</b> <b>E05</b> <b>E08</b>	HANDHELD AUTO DIAGNOSIS TEST FAILED	Contact technical service.
<b>E06</b>	GENERAL HANDHELD ERROR	Try to repeat exposure. If the problem persists, contact technical service.
<b>E07</b>	RF SIGNAL TOO LOW	Handheld lost the signal. Try to repeat exposure. If the problem persists, contact technical service.
<b>E09</b>	HANDHELD SERIAL NUMBER INCORRECT OR NOT INITIALIZED	Contact technical service.
<b>E10</b> <b>E12</b> <b>E13</b> <b>E16</b>	X-RAY UNIT INTERNAL ERROR	Contact technical service.

<b>E11</b>	COLLIMATOR SELECTION NOT CONSISTENT	After turning the rectangular collimator on or off, wait a few seconds to allow the icon on the handheld to be updated.
<b>E14</b> <b>E15</b>	GENERIC GENERATOR FAULT	Contact technical service.
<b>E17</b>	DEVICE OVERHEATED HEAD RELEASED	WAIT APPROXIMATELY 15 MINUTES FOR AUTOMATIC SYSTEM RESET
<b>E18</b> <b>E19</b>	SUPPLY VOLTAGE TOO HIGH/LOW	CHECK THE SUPPLY SYSTEM. IF THE PROBLEM PERSISTS, CALL TECHNICAL SUPPORT.
<b>E30</b>	INTERNAL ADJUSTMENT PROBLEM	REPEAT THE X-RAY. IF THE PROBLEM PERSISTS, CALL TECHNICAL SUPPORT.
<b>E31</b> <b>E32</b>	REMOTE CONTROL ERROR	REDUCE THE DISTANCE BETWEEN THE REMOTE CONTROL AND X-RAY HEAD AND THEN REPEAT THE X-RAY. FOLLOW THE INFORMATION GIVEN ON HOW TO PROPERLY USE THE HAND HELD'S ANTENNA.
<b>E33</b>	NO POWER SUPPLIED TO X-RAY GENERATOR	X-ray generator or arm cord may be faulty. Contact technical support service..

As regards the other error codes, **CONTACT** the technical service department.

## 7. Routine maintenance

### ATTENTION!

Any technical maintenance work required must be carried out by qualified personnel or by a specialised technician authorised by the manufacturer.

It is the user's responsibility to check that routine maintenance is carried out by an authorised technician at least every 2 years. The maintenance methods are specified in the Technical Service Manual possessed by the Authorised Technicians.

## 8. Cleaning and disinfection

The x-ray unit may be a source of infection passed from one patient to another if it is not cleaned properly. For this reason it should be disinfected on the outside every day after use.

If digital X-ray sensors are used make sure they are always used with disposable hygienic covers. When disinfecting the x-ray unit use soft disposable paper and avoid using corrosive substances or immersing the unit in liquids.

To avoid damaging the plastic materials use products containing:

- **Ethanol 96%:**  
Concentration: maximum 30 g per 100 g of disinfectant.
- **Propanole:**  
Concentration: maximum 20 g per 100 g of disinfectant.
- **Combination of ethanol and propanole:**  
Concentration: the combination of the two should be maximum 40 g per 100 g of disinfectant.

Compatibility tests between plastics and the following products have been carried out with no negative consequences:

- Incidin Spezial (Henkel Ecolab).
- Omnidid (Omnident)
- Plastisept (ALPRO) (not tuberculocide as not an alcohol-based disinfectant)
- RelyOn Virkosept (DuPont).
- Green & Clean SK (Metasys) (not tuberculocide as not an alcohol-based disinfectant)

### ATTENTION!

- **Do not use products containing isopropyl alcohol (2-propanol, iso-propanol).**
- **Do not use products that contain sodium hypochlorite (bleach).**
- Do not use cleaners that contain phenol.
- **Do not spray the selected products directly on the surfaces.**
- **Never combine products with each other or with liquids other than the products listed above.**
- **All products must be used as directed by the manufacturer.**

### Cleaning and disinfecting instructions.

Use soft disposable paper towels or sterile gauze to clean and disinfect. Do not use sponges or in any case, any material that can be reused.

### ATTENTION!

- **Shut off the unit prior to clean and disinfecting the external parts.**
- **Never lubricate the pivot point of the x-ray cone as proper operation of the locking system may be compromised.**
- **All material used to clean and disinfect must be thrown away.**

## 9. Dismantling the equipment when no longer used

As set out in Directives 2002/95/ EC, 2002/96/ EC and 2003/108/ EC, on the restrictions of the use of certain hazardous substances in electrical and electronic equipment along with collection, treatment, recycling and disposal of waste electrical and electronic equipment the latter must be treated as municipal waste, therefore sorted and collected separately. When new equipment of equivalent type is purchased the waste equipment should be returned to the distributor on a one-to-one basis for disposal.

As far as reuse, recycling and other forms of waste recovery mentioned above are concerned, the manufacturer is responsible for the actions specified by individual local laws.

Efficient collection of sorted waste separately to recycle and treat waste electrical and electronic equipment aids in preventing negative environmental impacts while protecting human health. In addition it facilitates recycling of the materials used to construct the equipment.

The crossed out wheeled bin placed on the equipment indicates that the waste equipment must be collected separately from other waste.


### ATTENTION!

Illegal waste disposal carries heavy fines defined by local laws.

## 10. Technical specifications

### SPECIFICATIONS

- Rated voltage: 230Vac/115Vac (according to the model).
- Max. mains voltage fluctuation:  $\pm 10\%$ .
- Maximum current: 6A for the 230Vac version; 10A for the 115Vac version at 60KV 7mA
- Frequency: 50/60Hz.
- Maximum absorbed power: 1,4KVA.
- Apparent line resistance: 0,5 $\Omega$ (230Vac), 0,2 $\Omega$ (115Vac).
- Fuses: 6.3A T for the 230Vac version; 10A T for the 115Vac version
- Generator: constant potential type.
- High nominal voltage: 60 / 63 / 65KV.
- Nominal current: 6 / 7mA.
- Power requirements at 0.1 sec: 420W (60KV 7mA), 378W (63KV 6mA), 390W (65KV 6mA).
- Current reference time: 0.7 mAs (7mA – 0.1s) / 0.6 mAs (6mA – 0.1s).
- Focal spot : 0,4mm
- Total filtration : 2.5mm Al @ 65KV.
- Half-value layer (HVL): >2mm Al @ 65KV.
- Leaked radiation: <0,25mGy / h at 1 metre from focusing point at 65KV 6mA, duty cycle 1:60
- Ability to be reproduced: 0,05.
- Electrical classification: Class I - Type B, intermittent service.
- Set exposure time: from 0.010 to 1,000 seconds
- Accuracy of times indicated:  $\pm 10\%$ .
- mGy display precision:  $\pm 30\%$ .


### WEIGHTS


- Weight of the unit with packaging: 38Kg max.
- Weight of the hand-held control panel: 0,3 Kg.

### RADIOGENIC TUBE F.s.=0,4 mm

- Radiogenic tube: TOSHIBA D-041.
- Focal spot: 0,4 mm in compliance with IEC 336 / 1993.
- Tolerance for position of the focal spot along the reference axis:  $\pm 2\%$ .
- Nominal high voltage and maximum allowable current: (65KV, 7mA)  $\pm 10\%$ .
- Anode construction material: Tungsten (W).
- Anode inclination: 12,5°.
- Anode thermal load: 4,3 KJ (6 KHU).
- Maximum continuous heat dissipation: 100 W.
- Operating cycle: 1:60 (1 second exposure - 60 seconds pause time).


### Anode Thermal Characteristics

## CONE TECHNICAL SPECIFICATIONS

- SSD=30cm (12"), X-ray beam less than or equal to Ø55mm.

A – REFERENCE AXIS

B – FOCAL SPOT IDENTIFICATION

## HANDHELD BATTERIES

- Type 2 x AA Alkaline 1.5 V

## TECHNICAL FACTOR MEASURE


The high voltage value is measured with a non-invasive instrument.



The anode current is controlled inside with measurement resistors and circuits to obtain very precise measurements. Operation of the circuits is checked at the time of testing. Once assembled, the anode current can no longer be directly measured.

The exposure time should be evaluated by measuring the time that elapses from the moment in which high voltage exceeds 75% of the nominal value to the moment in which it drops below this value. Considering the high gradient of the rising and trailing edges of the anode voltage and squaring due to inherent filtration, use of a threshold ranging from 25% to 75% may be considered non-influential.

## 11. Dimensions

All dimensions are expressed in millimeters (inches).

## 12. Identification nameplates

### WARNING!


Never remove the identification nameplates provided on the generator, central control unit and collimator cone.

#### Central control unit (NAMEPLATE)

The nameplate is found beside the main switch.

Details shown on the nameplate:

- Name of the manufacturer
- Name of the equipment
- Rated voltage
- Type of current
- Rated frequency
- Maximum absorbed power
- Serial number
- Manufacturing date


#### X-ray unit head

The nameplate is found on the lower cover at the back of the radiogenic unit.

Details shown on the nameplate:

- Name of the manufacturer
- Name of the equipment
- Specifications
- Model and serial number of x-ray tube
- Equipment serial number
- Date of manufacture


#### Collimator

The nameplate for the collimator is found outside it.

Details shown on the nameplate:

- Name of the manufacturer
- Type of cone
- Serial number
- Date of manufacture


#### Handheld

The nameplate for the handheld is found in the battery compartment.

Details shown on the nameplate:

- Name of the manufacturer
- Name of the equipment
- Rated voltage
- Number and type of batteries
- Serial number

**DICHIARAZIONE “CE” DI CONFORMITÀ  
“CE” CONFORMITY DECLARATION  
DECLARATION DE CONFORMITÉ “CE”  
ERKLÄRUNG VON “CE” ZUSTIMMUNG  
DECLARACION DE CONFORMIDAD “CE”  
DECLARAÇÃO DE CONFORMIDADE “CE”  
ΔΗΛΩΣΗ ΠΙΣΤΟΤΗΤΑΣ “CE”  
ДЕКЛАРАЦИЯ СООТВЕТСТВИЯ “CE”  
DEKLARACJA ZGODNOSCI WE  
“CE” UYGUNLUK BELGESİ**

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Imola, li \_\_\_\_\_

**Bussolari Paolo**  
**Managing Director**

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