

Veterinary Refractometer Handrefraktometer

726000

What do the ND, UG and S.P. refer to?

ND stands for refractive index scale. This scale goes from 1333–1360.

UG is the urine specific gravity scale and goes from 1000–1050 (be careful not to misinterpret it as going from 0–50).

SP is the serum protein scale and goes from 0–12 g/100 mL. Most normal values are calculated as g/L, therefore the value on the scale of the refractometer has to be multiplied by 10 to get g/L.

Was bedeuten ND, UG und S.P.?

ND steht für refractive index (Brechungsindex). Die Skala reicht von 1333–1360.

UG bezeichnet das Uringewicht. Die Skala reicht von 1000–1050 (Achtung: Die Skala kann leicht als 0–50 missinterpretiert werden).

S.P. ist das Serum-Protein. Die Skala reicht von 0–12 g/100 mL. Die meisten Normalwerte werden in g/L angegeben. Hierfür muss der abgelesene Wert mit dem Faktor 10 multipliziert werden.



EICKEMEYER®

GERMANY | Eltstraße 8 | 78532 Tuttlingen | Deutschland | T +49 7461 96 580 0 | F +49 7461 96 580 90 | info@eickemeyer.de | www.eickemeyer.de

SWITZERLAND | Sandgrube 29 | 9050 Appenzell | Schweiz | T 071 788 23 13 | F 071 788 23 14 | info@eickemeyer.ch | www.eickemeyer.ch

UNITED KINGDOM | 3 Windmill Business Village | Brooklands Close | Sunbury-on-Thames | Surrey | TW16 7DY | UK | T 020 8891 2007 | www.eickemeyer.co.uk

DENMARK | Solbakken 26 | Hammelev | 6500 Vojens | Danmark | T +45 70 20 50 19 | info@eickemeyer.dk | www.eickemeyer.dk

NETHERLANDS | Bellweg 44 | 4104 BJ Culemborg | Nederland | T +31 345 58 9400 | info@eickemeyer.nl | www.eickemeyer.nl

ITALY | Via G. Verdi 8 | 65015 Montesilvano (PE) | Italia | T 085 935 4078 | F 085 935 9471 | info@eickemeyer.it | www.eickemeyer.it

POLAND | Al. Jana Pawła II 27 | 00-867 Warszawa | Polska | T 22 185 55 76 | F 22 185 59 40 | info@eickemeyer.pl | www.eickemeyer.pl

CANADA | 617 Douro Street, Suite #205 | Stratford | Ont. | Canada | N5A 0B5 | T +1 519 273 5558 | F +1 519 271 7114 | info@eickemeyervet.ca | www.eickemeyercanada.ca