



Technical Specifications



Generator	32kW / max. 400mA
Voltage range	40 - 150kV
Tube	Focal spots 1,0 / 2,0 mm
Exposure time range	1 ms - 10000 ms
Programs	352 pre-programmed
Combi server	1 TB
Monitor	X-HW More Multitouch
Table top	160 cm x 81 cm, floating Longitudinal ±30 cm (11,8") Transversal ±10 cm (3,9 ")
Mains connection	230 V - 50 / 60 Hz 1-phase (wall plug)
Detector	X-DR XL wired, 43 cm x 43 cm (17" x 17"), CsI
X-AQS PACS Vet	X-AQS PACS Vet X-AQS Module Archiving on CD/DVD X-AQS Module Data Compression contains for Client X-AQS Module Patient CD/DVD X-AQS Module Import Bitmap Images X-AQS Module Import DICOM Images X-AQS Module Import AQS Images X-AQS Module Image Export X-AQS Module Mail Export
X-AQS Client	1 workstation (expandable up to 4 workstations) X-AQS Viewer X-AQS Diagnosis License per Client, floating license for X-AQS Software